

# APPLICATION

## Using a Step by Step Dilution to Improve $\beta$ -Gone<sup>TM</sup> $\beta$ -Glucuronidase Removal Methods

Matt Brusius and Jennifer Watson  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



**Matt Brusius**  
Product Manager  
*Matt Brusius is an avid ice hockey player. He likes skating backwards and taking slapshots from the point.*



### Overview

The standard protocol for  $\beta$ -Gone  $\beta$ -Glucuronidase Removal Products suggests diluting the sample to a final concentration of 40 % methanol prior to filtration. In this technical note we will explore how using a smaller dilution with a smaller percentage of methanol in a step by step approach, as opposed to diluting the solution all at one time, can yield acceptable results for most comprehensive drug research panel suites.

### Materials and Methods

All reagents and solvents were HPLC or analytical grade. Analyses were performed using an API 4000<sup>TM</sup> LC/MS/MS (SCIEX, Framingham, MA)

### Sample Preparation

#### Enzymatic Hydrolysis

IMCSzyme<sup>®</sup> Hydrolysate Mix was prepared as follows:

1. Add 10  $\mu$ L of analyte spike (1  $\mu$ g/mg) to 140  $\mu$ L of urine
2. Dilute with 80  $\mu$ L of IMCS buffer (IMCS Part No.: 04-EZRHB-20) and add 30  $\mu$ L IMCS Enzyme (IMCS Part No.: 04-EIF-010) then vortex for 15 seconds
3. Proceed to  $\beta$ -Glucuronidase Removal Protocols

#### $\beta$ -Glucuronidase Removal

##### $\beta$ -Gone Recombinant Enzyme – 20 % Methanol Step by Step Addition

1. Dilute 200  $\mu$ L urine hydrolysate and load onto the Recombinant  $\beta$ -Gone 96-Well Plate (Phenomenex Part No.: 8E-S139-TGA) then apply 2-5" Hg vacuum for one minute
2. Add 67  $\mu$ L of 1 % formic acid in methanol and apply 5" Hg vacuum for one minute.
3. Add 67  $\mu$ L of 1 % formic acid in water and apply 5" Hg vacuum for one minute
4. Collect pooled eluent and inject 10  $\mu$ L

### Results and Discussion

**Table 1** provides absolute recovery values of the samples prepared as described in the "Materials and Methods" section. The undiluted urine hydrolysate is loaded directly onto the well plate and diluted by 67  $\mu$ L of 0.1% formic acid in methanol. The dilution is then repeated with 67  $\mu$ L of 0.1% formic acid in water. The order of the dilution is critical and the methanolic mixture must be added first to effectively disrupt hydrophobic interaction which improves recovery and reproducibility for more non-polar species (i.e. THC-COOH). The secondary elution of 0.1% formic acid in water is required to ensure that the entire bed is sufficiently saturated. 133  $\mu$ L is equal to approximately 4 bed volumes of the 30 mg sorbent bedmass, which is the established practical minimum volume for wash and elution conditions in traditional Solid Phase Extraction (SPE).


Data not shown in this experiment also suggests that the same protocol can be used for the Non-Recombinant version of  $\beta$ -Gone, but requires a secondary elution of 117  $\mu$ L of 0.1% formic acid in water to fully saturate the larger bed mass. This results in a 33% absolute recovery of THC-COOH, but with a CV of 9% (n=8).



**Table 1.**  
Recovery for Recombinant  $\beta$ -Gone<sup>TM</sup> using a step by step dilution

Analyte	Average % Recovery	%CV (N=8)
6-MAM	105	4
Alprazolam	97	6
Buprenorphine	76	13
Codeine	91	8
Diazepam	88	11
Fentanyl	88	7
Flurazepam 1	96	6
Hydrocodone	83	8
Hydromorphone	107	16
MDMA	89	8
Methadone	95	9
Methamphetamine	105	5
Morphine	83	17
Naloxone	105	10
Norbuprenorphine	106	10
Nordiazepam	101	7
Norfentanyl	106	9
Oxycodone	97	11
PCP	91	11
Temazepam	88	9
THC-COOH	55	8

**Figure 1:**  
Chromatogram of  $\beta$ -Gone with step by step dilution



#### HPLC Conditions

**Column:** Kinetex<sup>®</sup> 2.6  $\mu$ m Biphenyl MercuryMS<sup>™</sup> Cartridge  
**Dimensions:** 20 x 2.0 mm  
**Part No.:** 00M-4622-B0-CE  
**Mobile Phase:** A: 0.1 % Formic acid in Water  
 B: 0.1 % Formic acid in Methanol  
**Gradient:**

Time (min)	B (%)
0	20
2	100
2.5	100
2.51	20
4	20

**Flow Rate:** 0.7 mL/min  
**Detection:** MS/MS (API 4000<sup>™</sup>, SCIEX)

**Ordering Information**

**β-Gone™ β-Glucuronidase Removal Products**

Part No.	Description	Unit
8B-S139-TAK	1 mL Tubes, Recombinant Enzyme	100/Box
8B-S322-DAK	1 mL Tubes, Non-Recombinant Enzyme	100/Box
8E-S139-TGA	96-Well Plate, Recombinant Enzyme	1/Box
8E-S322-DGA	96-Well Plate, Non-Recombinant Enzyme	1/Box
8E-S322-TUK	2 mL Centrifuge Tubes, Recombinant and Non-Recombinant Enzyme	100/Box



LC/MS Cartridges		
2.6 μm	Phase	20 x 2.0 mm
Kinetex	Biphenyl	00M-4622-B0-CE

**Standard Cartridge Holders**

Part No.	Description
CH0-5846	10 mm standard holder
CH0-5845	20 mm standard holder

**Vacuum Manifolds**

Part No.	Description	Unit
<b>12-Position Vacuum Manifold for Tubes*</b>		
AH0-6023	12-Position Vacuum Manifold Set, complete assembly	ea
<b>24-Position Vacuum Manifold for Tubes*</b>		
AH0-6024	24-Position Vacuum Manifold Set, complete assembly	ea
<b>96-Well Plate Manifold</b>		
AH0-8950	96-Well Plate Manifold, Universal with vacuum gauge	ea



\* Manifolds include: Vacuum-tight glass chamber, vacuum gauge assembly, polypropylene lid with gasket, male and female luers and yellow end plugs, stopcock valves, collection rack assemblies, polypropylene needles, lid support legs. Waste container included with 12-position manifold.

**Presston™ 100 Positive Pressure Manifold**

Part No.	Description
AH0-9334	Presston 100 Positive Pressure Manifold, 96-Well Plate
AH0-9342	Presston 100 Positive Pressure Manifold, 1 mL Tube Complete Assembly
AH0-9347	Presston 100 Positive Pressure Manifold, 3 mL Tube Complete Assembly
AH0-9343	Presston 100 Positive Pressure Manifold, 6 mL Tube Complete Assembly



The Presston 100 96-Well Positive Pressure Manifold can also process 1, 3, and 6 mL tubes using the following adapter kits

**Presston 100 Tube Adapter Kits (for AH0-9334)**

Part No.	Description
AH0-9344	1 mL Tube Adapter Kit
AH0-9345	3 mL Tube Adapter Kit
AH0-9346	6 mL Tube Adapter Kit



WARRANTY Phenomenex warrants that for a period of 12 months following delivery, the Presston 100 Positive Pressure Manifold you have purchased will perform in accordance with the published specifications and will be free from defects in materials or workmanship. In the event that the Presston 100 Positive Pressure Manifold does not meet this warranty, Phenomenex will repair or replace defective parts. Please visit [www.phenomenex.com/Presston](http://www.phenomenex.com/Presston) for complete warranty information.



# APPLICATION

## Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auiinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

## China

t: +86 (0)20 2282-6668  
f: +86 (0)20 2809-8130  
chinainfo@phenomenex.com

## Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

## India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

## Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

## All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and demensions, return the products with comparative data within 45 days for a FULL REFUND.

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

### Trademarks

Kinetex is a registered trademark and  $\beta$ -Gone, MercuryMS, and Presston are trademarks of Phenomenex. IMCSymze is a registered trademark of Integrated Micro-Chromatography Systems, LLC. API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

© 2016 Phenomenex, Inc. All rights reserved.