








pH HPLC METHOD DEVELOPMENT



Adjust pH to Alter Selectivity



Column: Gemini-NX 5 μ m C18
 Dimensions: 150 x 4.6 mm
 Part No.: 00F-4454-E0
 Mobile Phase: (see chromatograms)
 Gradient: A/B (95:5) to (5:95) for 10 min, hold for 2 min
 Flow Rate: 1.5 mL/min
 Temperature: Ambient
 Detection: UV @ 254 nm
 Sample
 Analytes:
 1. Amitriptyline
 2. Naproxen
 3. Toluene



- ### Method Development Tips
- Compounds being analyzed should be neutral for reversed phase chromatography
 - For acidic and basic compounds, keep the mobile phase 2 pH units above or below the pK_a value
 - For increased column lifetime, use a < 20 mM buffer concentration
 - Be sure to stay within the buffering range of the buffer being used



Retention of Different Functional Groups

Analytes contain functionalities which will affect retention on a reversed phase column. Whether a functionality is protonated or deprotonated will usually affect the chromatographic behavior of a compound; altering the pH of a separation will often change an elution profile. A general observation is that an acidic functionality will tend to have greater retention and efficiency at a pH below its pK_a value and less retention at a pH above its pK_a value. Basic compounds follow the opposite trend: basic functionalities are often less retained at pH's below their pK_a values and demonstrate greater retention and better peak shape at pH's above their pK_a values.

Protonated Functional Group	Deprotonated Functional Group	Approx. pK_a
<chem>OP(=O)(O)O</chem>	<chem>OP(=O)(O)[O-]</chem>	1.5
<chem>OS(=O)(=O)O</chem>	<chem>OS(=O)(=O)[O-]</chem>	2.0
<chem>OC(=O)O</chem>	<chem>OC(=O)[O-]</chem>	4.7
<chem>C1=CC=NC=C1</chem>	<chem>C1=CC=NC=C1</chem>	5.3
<chem>C1=CN=CN=C1</chem>	<chem>C1=CN=CN=C1</chem>	7.0
<chem>CN(C)C</chem>	<chem>CN(C)C</chem>	9.8
<chem>Oc1ccccc1</chem>	<chem>[O-]c1ccccc1</chem>	10.0
<chem>C[NH3+]</chem>	<chem>CN</chem>	10.6
<chem>C1CCNCC1</chem>	<chem>C1CCNCC1</chem>	10.7

Buffer Selection For Method Development

Next to organic solvent, buffer selection is the most important variable in HPLC method development. There are several criteria in selecting the appropriate buffer. The first is choosing a buffer that has a pK_a near the desired working pH. Other criteria like ionic strength, ion-pairing properties, as well as mass spectrometer compatibility should be considered before selecting any mobile phase modifier.

Buffer†	pK_a	Buffer Range (pH)	MS Compatible
Trifluoroacetic Acid	< 2	< 2.5	•‡
Phosphoric Acid (pK_1)	2.1	1.1 - 3.1	
Citric Acid (pK_1)	3.1	2.1 - 4.1	
Formic Acid	3.8	2.8 - 4.8	•
Citrate (pK_2)	4.7	3.7 - 5.7	
Acetic Acid	4.8	3.8 - 5.8	•
Citrate (pK_3)	5.4	4.4 - 6.4	
Carbonate (pK_1)	6.4	5.4 - 7.4	•
Phosphate (pK_2)	7.2	6.2 - 8.2	•
Triethanolamine	7.8	6.8 - 8.8	•
TRIS	8.3	7.3 - 9.3	
Diethanolamine	8.9	7.9 - 9.9	•
Ammonia	9.2	8.2 - 10.2	•
Ethanolamine	9.5	8.5 - 10.5	•
Carbonate (pK_2)	10.3	9.3 - 11.3	•
Diethylamine	10.5	9.5 - 11.5	•
Triethylamine	11.0	10.0 - 12.0	•
Piperidine	11.1	10.1 - 12.1	
Phosphate (pK_3)	12.3	11.3 - 13.3	

† Common buffers are in bold.
 ‡ TFA can be used at low concentrations for LC/MS applications but can affect MS sensitivity.



- Gemini[®] pH-LC[™]**
 - pH 1-12 Stability
 - High Efficiency
 - Rugged and Reproducible Results
- Gemini-NX**
 - Longer Lifetime Under Extreme pH Conditions

www.phenomenex.com/gemini
 Phenomenex products are available worldwide. Email us at info@phenomenex.com.

