

APPLICATIONS

Effect of Alkylamine Based Ion-pairing Reagents in Analysis of siRNA

Ivan Lebedev, Brian Rivera, and Namrata Saxena
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Silencing RNA, or siRNA, are short, double-stranded RNA molecules which are not only an important research tool in molecular biology, but also an emerging therapeutic modality. Numerous chromatographic methods have been developed for analysis of these double stranded RNA molecules, including LC-ESI-MS and LC-UV methods.

Initial oligonucleotide LC-MS methods used high mounts of TEA and HFIP (16 and 400 mM, respectively) to facilitate the retention of these negatively charged, polar compounds¹. However, with modern HPLC media, this mobile phase composition may not be optimal for both chromatographic performance and electrospray ionization. Although many alkylamine ion-pair have been investigated for MS optimization², in this application note, we show that chromatographic separation of siRNA sense and anti-sense strands may be more optimal with TEA.

Figure 1 shows a 21mer double-stranded RNA using TEA as the ion-pairing reagent. Even with 12.5 mM HFIP, the relatively low concentration 4 mM TEA provides sufficient retention of the siRNA, with good resolution of sense and anti-sense strand; although not baseline, this would be acceptable to obtain high quality spectral data if characterizing each strand separately by high resolution MS.

N,N-diisopropylethylamine (DIEA) is another commonly used alkylamine, which often provides superior chromatographic separation and improved electrospray desorption for oligos. However, as seen in **Figure 2** for the same 21mer double-stranded RNA, we see some aberrant chromatography. Although there are seemingly more earlier eluting impurities, loss of resolution sense and anti-sense is observed. Further MS characterization would need to be performed to determine the actual impact on analysis. Similar chromatographic behavior is observed with more hydrophobic alkylamine, including N,N-dimethylcyclohexylamine (DMCHA) and hexylamine (HA). After optimizing ion pair concentration and gradient profile, separation using DMCHA is not favorable (**Figure 3**), and this odd chromatographic behavior is even more pronounced with HA (**Figure 4**).

It is important to note that if the intent is MS sensitivity, TEA may not perform as well as other alkylamine, as reported in the literature. Further, the more hydrophobic alkylamine show more impurities, though without the MS characterization it is unclear if these are indeed impurities/failure sequences or if these impurity peaks are a result of chromatographic aberration. Replicate injections showed consistent chromatography, though day to day data showed different chromatography (data not shown).

In summary, if the separation of sense and anti-sense strand when characterizing double-stranded oligonucleotides is prioritized, TEA will likely yield more favorable chromatography and more reproducible results. As such, other factors impacting MS performance, including modulation of HFIP concentration, might need to be implemented.

LC-UV Conditions

Column: bioZen™ 2.6 µm Oligo
Dimension: 100 x 2.1 mm
Part No.: [00D-4790-AN](#)
Mobile Phase: A: 12.5 mM HFIP, IP Reagent as mentioned in Water
B: 12.5 mM HFIP, IP Reagent as mentioned in MeOH
Gradient: As mentioned
Flow Rate: 0.3 mL/min
Injection: 2 µL (1 µg/mL)
Temperature: 55 °C
Detection: UV-Vis @ 260 nm
Sample: siRNA

¹Gilar M, Fountain KJ, Budman Y, et al. Ion-pair reversed-phase high-performance liquid chromatography analysis of oligonucleotides: retention prediction. *J Chromatogr A*. 2002;958(1-2):167-182

²Basiri, Babak et al. "Assessing the Interplay between the Physicochemical Parameters of Ion-Pairing Reagents and the Analyte Sequence on the Electrospray Desorption Process for Oligonucleotides." *Journal of the American Society for Mass Spectrometry* vol. 28,8 (2017): 1647-1656.

Figure 1. Separation of 21mer double-stranded RNA, 4 mM TEA




Figure 2. 21mer double-stranded RNA, 2 mM DIEA




Figure 3. 21mer double-stranded RNA, 0.5 mM DMCHA





Figure 4. 21mer double-stranded RNA, 0.5 mM HA



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com