

APPLICATION

Chiral Separation of the Drug Product Levetiracetam Using the United States Pharmacopeia Method on the Latest L51 USP Column - Lux[®] Amylose-1

Emily Gish, Michael McCoy, Marc Jacob and Sean Orlowicz
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Emily Gish
Research Associate
Emily is an avid baker. She loves throwing themed parties to entertain friends.

In this technical note, we report the enantiomeric separation of Levetiracetam using the approved United States Pharmacopeia method. We perform the chiral analysis on the latest L51 USP column, the polysaccharide-based chiral stationary phase Lux Amylose-1 column, and compare analytical result as well as performance to CHIRALPAK[®] AD-H[®] column.

Introduction

The United States Pharmacopeia (USP) establishes methods and provides reference standards for medicines, food ingredients, dietary supplement products, and active pharmaceutical ingredients (API). These USP methods and standards are used by regulatory agencies and pharmaceutical manufacturers to ensure that APIs or drug products are of the appropriate identity, as well as strength, quality, purity, and consistency. In the case of chiral drugs, such as Levetiracetam (**Figure 1**.) a chiral HPLC method is generally used to assess chiral purity.

Levetiracetam (trade name Keppra) is an anticonvulsant medication mainly used to treat seizure disorder (epilepsy), as well as psychiatric and neurologic conditions in humans and animals. As of November 2008, the drug is available as a generic brand in the United States and the United Kingdom.

Levetiracetam is the biologically active S-enantiomer of the racemic compound Etiracetam. A USP method using a L51 USP column is published to perform enantiomeric separation of a racemic mixture and quantify the R-enantiomer in the drug product Levetiracetam.

In this technote, we report the results for the limit of Levetiracetam (R-enantiomer) following the published USP method using the newest L51 USP column, Lux Amylose-1, and compare its performance to the well known CHIRALPAK AD-H (L51) column.

Materials and Methods

All analyses were performed using a Shimadzu SCL-10AVP LC system (Shimadzu Corporation, Canby, OR, USA) equipped with binary pump, in-line degasser, multi-wavelength UV detector and autosampler. The Lux Amylose-1 column used for analysis was obtained from Phenomenex (Torrance, CA, USA) and the CHIRALPAK AD-H column used was obtained from DAICEL[®] Corporation (Fort Lee, NJ, USA). USP Reference Standards were purchased from Sigma-Aldrich (St. Louis, MO, USA). All solvents were purchased from Honeywell (Morristown, NJ, USA).

USP sample and definition

System suitability solution: 0.1 mg/mL of USP Levetiracetam (RS racemic mixture) in mobile phase


Sample solution: 10mg/mL of USP Levetiracetam (S-enantiomer) in mobile phase

Each sample solution was subjected to the experimental conditions described below, first using Lux Amylose-1, followed by CHIRALPAK AD-H.

HPLC Conditions

- Column:** Lux 5 μ m Amylose-1
CHIRALPAK 5 μ m AD-H
- Dimensions:** 250 x 4.6 mm
- Mobile Phase:** Hexane/Ethanol (80:20)
- Flow Rate:** 1 mL/min
- Temperature:** 25 $^{\circ}$ C
- Detection:** UV @ 215 nm
- Sample:** 1. Levetiracetam R-enantiomer
2. Levetiracetam S-enantiomer (biologically active compound)

Figure 1.
Chemical Structure of Levetiracetam



Results and Discussion

The USP HPLC analytical method for the limit of Levetiracetam (R-enantiomer) was evaluated on both the Lux[®] Amylose-1 and CHIRALPAK[®] AD-H[®] columns to compare product performance. The USP standards were dissolved in the mobile phase and injected on a Shimadzu HPLC system as described in the material and methods section.

The USP method gives guidelines for the relative retention time (RRT) of the Levetiracetam enantiomers in the system suitability section. According to the monograph, the R-enantiomer and S-enantiomer RRT should be 0.55 and 1.0 respectively for the system suitability solution.


The resolution (R_s) criteria for the system suitability solution between the two enantiomers of Levetiracetam must be no less than (NLT) 4.0 and the amount of R-enantiomer shall be no more than (NMT) 0.8 % for the drug to meet the chiral purity acceptance criteria.

USP reference standards were first injected on the Lux Amylose-1 column following the USP monograph method. The retention times and resolution of the system suitability solution using the Lux Amylose-1 column are shown in **Figure 2a**. The retention times for the R- and S-enantiomers were 8.72 minutes (RRT 0.52) and 16.78 (RRT 1.0) minutes, respectively. The resolution obtained with the Lux Amylose-1 column was 13.9, well within the USP method criteria of no less than 4.0. Retention times and resolution for system suitability solution on CHIRALPAK AD-H column are shown in **Figure 2b**. The retention times for the R- and S-enantiomers were 8.66 minutes (RRT 0.52) and 16.61 (RRT 1.0) minutes, respectively. The results using the CHIRALPAK AD-H column are similar to Lux Amylose-1 for the RRT results but resolution was slightly lower with 13.4 vs 13.9 for the Lux column. The higher resolution is related to the better column efficiency of the Lux Amylose-1.

Figure 3a displays the chromatogram of Levetiracetam sample solution run on Lux Amylose-1 column. The calculated percentage of R-enantiomer is 0.29 %, well below the 0.8 % required by the USP method. **Figure 3b** shows the chromatogram of Levetiracetam sample solution run on CHIRALPAK AD-H column. The calculated percentage of R-enantiomer is 0.58 %, also adhering to the criteria of NMT 0.8 %. **Table 1** displays the corresponding results on both the Lux Amylose-1 and the CHIRALPAK AD-H columns for percentage of R-enantiomer in the sample solution.


Figure 2. Levetiracetam system suitability on Lux Amylose-1 and CHIRALPAK AD-H

Figure 2a. System suitability solution on Lux Amylose-1



Peak	Ret Time	Area	Plates/m	R_s	Area%
1	8.72	765671	37686	-	48.42
2	16.78	815578	61626	13.90	51.58

Figure 2b. System suitability solution on CHIRALPAK AD-H




Peak	Ret Time	Area	Plates/m	R_s	Area%
1	8.66	780313	34961	-	50.15
2	16.61	775561	58599	13.39	49.85

* Comparative separations may not be representative of all applications


Figure 3.
Levetiracetam sample solution on Lux[®] Amylose-1 and CHIRALPAK[®] AD-H[®]

Figure 3a.
Sample solution on Lux Amylose-1



Peak	Ret Time	Area	Plates/m	Rs	Area%
1	8.67	448388	35398	-	0.23
2	15.44	194951304	6363	5.47	99.77

Figure 3b.
Sample solution on CHIRALPAK AD-H



Peak	Ret Time	Area	Plates/m	Rs	Area%
1	8.66	453801	35469	-	0.21
2	14.93	215663523	4568	4.6	99.80

Table 1.
Percentage of the R-enantiomer in Levetiracetam sample solution on Lux Amylose-1 and CHIRALPAK AD-H.

USP L51 Column	ru	rs	cu	cs	%R
Lux Amylose-1	448388	765671	0.05	10	0.29
CHIRALPAK AD-H	453801	780313	0.05	10	0.58

The percentage of unwanted R-enantiomer in Levetiracetam sample solution is calculated using the given formula:

$$\%R\text{-enantiomer} = [(ru/rs) \times (cs/cu)] \times 100$$

ru is the peak response of the R-enantiomer in the sample solution, rs is the peak response of the R-enantiomer in the standard solution, cs is the concentration of the standard solution, and cu is the concentration of the sample solution

Conclusion

The results shown in this technote demonstrate that the newest L51 USP column, Lux Amylose-1, produces equivalent results when compared to the CHIRALPAK AD-H column for the enantiomeric separation of Levetiracetam following the USP monograph for enantiomeric purity.

The resolution of the system suitability solution on both columns adhere to the criteria stated in the approved USP method for the Levetiracetam enantiomers. The resolution between the two enantiomers was calculated to be 13.9 on the Lux Amylose-1 column which is a better result than the CHIRALPAK AD-H column with a calculated resolution of 13.4. The difference is mainly explained by the better efficiency seen with the Lux Amylose-1 column.

Finally, the percentage of the R-enantiomer calculated from prepared standard and sample solutions was 0.29 % on the Lux Amylose-1 column and 0.58 % on the CHIRALPAK AD-H column, well below the USP criteria of NMT 0.8 %.

* Comparative separations may not be representative of all applications



APPLICATION

Lux[®] Ordering Information

5 µm Analytical Columns (mm)					SecurityGuard [™] Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
Amylose-1	00B-4732-E0	00D-4732-E0	00F-4732-E0	00G-4732-E0	AJ0-9336
					for ID: 3.2–8.0 mm

5 µm Semi-Prep Columns (mm)		SecurityGuard Cartridges (mm)
Phases	250 x 10.0	10 x 10.0 [‡]
Amylose-1	00G-4732-N0	AJ0-9344
		for ID: 9–16 mm

5 µm Axia [™] Packed Preparative Columns (mm)					SecurityGuard Cartridges (mm)	
Phases	150 x 21.2	250 x 21.2	250 x 30	250 x 50	15 x 21.2**	15 x 30.0*
Amylose-1	00F-4732-P0-AX	00G-4732-P0-AX	00G-4732-U0-AX	00G-4732-V0-AX	AJ0-9338	AJ0-9339
					for ID: 18–29 mm	30–49 mm

* SecurityGuard Analytical Cartridges require holder, Part No. : KJ0-4282

‡ SemiPrep SecurityGuard[™] Cartridges require holder, Part No.: AJ0-9281

**HPLC PREP SecurityGuard Cartridges require holder, Part No. : AJ0-8223

SFC PREP SecurityGuard Cartridges require holder, Part No. : AJ0-8617

* HPLC PREP SecurityGuard Cartridges require holder, Part No. : AJ0-8277

SFC PREP SecurityGuard Cartridges require holder, Part No. : AJ0-8618

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Lux is a registered trademark, Axia and SecurityGuard are trademarks of Phenomenex. DAICEL, CHIRALPAK, and AD-H are registered trademarks of DAICEL Corporation. All such trademarks are used by Chiral Technologies under license from DAICEL Corporation. Chiral Technologies, Inc. is a subsidiary of DAICEL Corporation.

Disclaimer

Comparative separations may not be representative of all applications. Columns used for comparison were manufactured by DAICEL Corporation. Phenomenex is in no way affiliated with DAICEL Corporation.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

© 2015 Phenomenex, Inc. All rights reserved.

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com