

# APPLICATIONS

## Fast Analysis of Alcohol-Based Sanitizers by GC-FID using Zebron<sup>TM</sup> ZB-WAXPLUS<sup>TM</sup> GC Columns

Ramkumar Dhandapani and Rola El abaji  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview


Hand sanitizer consumer products typically come in the form of liquid or gel and are thought to be generally effective at killing microorganisms and decreasing infectious agents on our hands. Alcohol-based sanitizers are most common, with versions that contain 60% to 95% alcohol. These versions typically contain some combination of isopropyl alcohol, ethanol (ethyl alcohol), and/or n-propanol.

Presented in **Figure 1** is the separation of alcohols commonly tested for in alcohol-based sanitizers including methanol, ethanol, propanol and butanol, on a 100% aqueous stable ZB-WAXPLUS GC column. In addition to the fast 15-minute separation, high injection-to-injection reproducibility can be noticed for this challenging aqueous matrix. This high-throughput separation example employed a high percentage ethanol to mimic alcohol-based sanitizer. **Figure 2** represents separation of trihydric alcohols, such as glycerols, that are added as moisturizer in some versions of hand sanitizers. A run time of less than 15 minutes is extremely helpful for faster testing and expedited batch release based on the GC test results.

### GC-FID Conditions


<b>Column:</b>	Zebron ZB-WAXPLUS
<b>Dimension:</b>	30 meter x 0.25 mm x 0.25 µm
<b>Part No.:</b>	<a href="#">7HG-G013-11</a>
<b>Injection:</b>	30:1 @ 140 °C, 0.2 µL
<b>Recommended Liner:</b>	Zebron PLUS Single Taper Z-Liner <sup>TM</sup>
<b>Liner Part No.:</b>	<a href="#">AG2-4B13-05</a> ; <a href="#">AG2-0A13-05</a>
<b>Carrier Gas:</b>	Helium @ 1.4 mL/min (constant flow)
<b>Oven Program:</b>	35°C for 5 min to 85°C @ 10°C/min to 200°C @ 25°C/min for 1 min
<b>Detector:</b>	Flame Ionization (FID) @ 200 °C
<b>Sample:</b>	1. Acetaldehyde 2. Ethyl acetate 3. Methanol 4. Ethanol 5. Propanol 6. Isobutanol 7. 2-Methylbutanol 8. 3-Methylbutanol

**Figure 1: Separation of Alcohols in Aqueous Matrix in Less than 15 Minutes on a Zebron ZB-WAXPLUS GC Column**



App ID 15817

**Figure 2: Separation of Glycerol in Less Than 11 Minutes on a Zebtron ZB-WAXPLUS GC Column**



#### GC-FID Conditions

- Column:** Zebtron ZB-WAXPLUS  
**Dimension:** 30 meter x 0.25 mm x 0.25  $\mu$ m  
**Part No.:** [7HG-G013-11](#)  
**Injection:** Split 20:1 @ 225  $^{\circ}$ C, 0.5  $\mu$ L  
**Recommended Liner:** Zebtron PLUS Single Taper Z-Liner  
**Liner Part No.:** [AG2-0A13-05](#)  
**Carrier Gas:** Helium @ 3 mL/min (constant flow)  
**Oven Program:** 80  $^{\circ}$ C to 180  $^{\circ}$ C @ 20  $^{\circ}$ C/min for 20 min  
**Detector:** Flame Ionization (FID) @ 250  $^{\circ}$ C  
**Sample:** 1. Benzene  
 2. Acetonitrile  
 3. 1,2-Propanediol  
 4. Ethylene glycol  
 5. 1,3-Propanediol (internal standard)  
 6. Diethylene glycol  
 7. Glycerol  
 8. Triethylene glycol

# APPLICATIONS

## Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

### Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Poland

t: +48 (12) 881 0121  
pl-info@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

### All other countries Corporate Office USA

t: +1 (310) 212-0555  
info@phenomenex.com

## [www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Zebron, WAXPLUS, and Z-Liner are trademarks of Phenomenex.  
FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.