

APPLICATIONS

Analysis and Quantitation of Pre-Workout Powder Components in Matrix by LC/UV using Luna[®] 3 μ m HILIC HPLC Column

Morgan Kramer, Matthew Trass, Allen Misa, and Scott Krepich
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Scott Krepich
Senior Field Application Scientist
Scott enjoys surfing and eating.
He is crazy about chromatography,
because his mom is really into CSI
and thinks that is what he does.



Introduction

With the implementation of the food safety modernization act (FSMA) and new label claim verification requirements, attention to pre-workout supplements have increased in recent years. This newfound popularity has led to an increase in interest for the identification and quantitation of the key ingredients promoted on the labels for these pre-workout drinks. The understanding of what is contained within a consumable supplement is important to ensure that the product being consumed is accurately represented for both product safety and quality. In this work, we present a fast and highly efficient method for the identification and quantitation of four key ingredients of pre-workout powders using HILIC mode HPLC-UV/Vis analytical techniques.


Materials and Methods

Analytical standards of Creatine, Creatinine, Caffeine and Dicyandiamide were purchased from Sigma-Aldrich, and were used for method development purposes. Each of these components were then made into stock solutions and injected at various concentrations to develop a standard calibration curve (**Figures 3-6**).

Sample Preparation

10 mg of Pre-Workout powder was dissolved into 1 mL of DI water and vortexed for 1 minute at 3000 rpm to ensure the sample was fully dissolved. The samples were then diluted to a concentration of 2.5 mg/mL with DI water for HPLC injection.

Figure 1.
Optimized Separation of 4 Key Components



HPLC Conditions


Column: Luna[®] 3 μ m HILIC 200 Å
Dimensions: 150 x 4.6 mm
Part No.: 00F-4449-E0
Mobile Phase: A: 1 mM Sodium Acetate pH 5.5
 (pH Adjusted with acetic acid)
 B: Acetonitrile

Gradient:	Time (min)	% B
	0	90
	4	71
	6	49

Flow Rate: 1.5 mL/min
Temperature: 35°C
Detector: UV/Vis @ 218 nm (Bandwidth 4nm) (ambient)
HPLC System: Agilent[®] 1100 HPLC with Quaternary Pump
Injection: 1 μ L (in water)
Analytes: 1. Caffeine (200 ppm)
 2. Dicyandiamide (200 ppm)
 3. Creatinine (400 ppm)
 4. Creatine (400 ppm)

APPLICATIONS

Figure 2.
Optimized method using variable standard concentrations in water



HPLC Conditions

Column: Luna[®] 3 μm HILIC 200 Å
Dimensions: 150 x 4.6 mm
Part No.: 00F-4449-E0
Mobile Phase: A: 1 mM Sodium Acetate pH 5.5 (pH Adjusted with acetic acid)
 B: Acetonitrile


Gradient: Time (min)	% B
0	90
4	71
6	49

Flow Rate: 1.5 mL/min
Temperature: 35 °C
Detector: UV/Vis @ 218 nm (Bandwidth 4 nm) (ambient)
HPLC System: Agilent[®] 1100 HPLC with Quaternary Pump
Injection: 1 μL (in water)
Analytes: 1. Caffeine (200 ppm)
 2. Dicyandiamide (200 ppm)
 3. Creatinine (400 ppm)
 4. Creatine (400 ppm)

App ID 23861

Blue: 400/200 ppm Red:200/100 ppm Green: 100/50 ppm Pink: 50/25 ppm Brown: 25/12.5 ppm
 *X/Y Refers to the concentration of X (Creatine and Creatinine) and the concentration of Y (Caffeine and Dicyandiamide) in the mixture

Figure 3.
Calibration Curve for Caffeine 12.5 to 200 ppm



APPLICATIONS

Figure 4.
Calibration Curve for Dicyandiamide 12.5 to 200 ppm




Figure 5.
Calibration Curve For Creatinine 25 to 400 ppm






Figure 6.
Calibration Curve For Creatine 25 to 400 ppm



APPLICATIONS

Figure 7.
C4 Pre-workout 1µL injection of 2.5 mg/mL



HPLC Conditions

Conditions for both separations:

Column: Luna[®] 3 µm HILIC 200 Å

Dimensions: 150 x 4.6 mm

Part No.: 00F-4449-E0

Mobile Phase: A: 1 mM Sodium Acetate pH 5.5
(pH Adjusted with acetic acid)
B: Acetonitrile

Gradient: Time (min)	% B
0	90
4	71
6	49

Flow Rate: 1.5 mL/min

Temperature: 35°C

Detector: UV/Vis @ 218 nm (Bandwidth 4 nm) (ambient)


HPLC System: Agilent[®] 1100 HPLC with Quaternary Pump

Injection: 1 µL (in water)

- Analytes:**
1. Caffeine (200 ppm)
 2. Dicyandiamide (200 ppm)
 3. Creatinine (400 ppm)
 4. Creatine (400 ppm)

App ID 23862

Figure 8.
N.O. Xplode Pre-workout 1µL injection of 2.5 mg/mL



App ID 23863

APPLICATIONS

Figure 9.
Bang Pre-workout 1µL injection of 2.5 mg/mL





Figure 10.
Outlift Pre-workout 1µL injection of 2.5 mg/mL



APPLICATIONS

Table 1. Concentrations of matrix components present in a 1 μ L injection of a 2.5 mg/mL matrix in water (0.0025 mg injection)

Concentrations of Components (ppm)				
Sample	Caffeine	Dicyandiamide	Creatinine	Creatine
C4	40.2	0.0	8.1	162.1
N.O. Xplode™	33.1	0.8	0.1	223.8
Bang®	36.9	0.8	0.9	331.4
Outlift™	21.7	0.8	3.4	175.3

Table 2. Scaled up concentrations of matrix components accounting for serving size

Mass concentration of matrix components (mg)					
Sample	Caffeine	Dicyandiamide	Creatinine	Creatine	Serving Size (mg)
C4	104.54	0.00	21.08	421.42	6500
N.O. Xplode	244.83	5.73	0.56	1656.40	18500
Bang	389.62	8.97	9.69	3499.43	26400
Outlift	222.54	7.71	34.60	1795.38	25600

Results and Discussion

With the recent increased interest in the analysis of Pre-Workout powders, a rapid and robust analytical technique was needed for the fast and efficient analysis of their key components. Due to the high polarity of the compounds present in these matrices, traditional reversed phase HPLC techniques can become rather difficult and time consuming, and in many cases, could even result in inaccurate quantitation and poor peak shape. Here we presented a highly selective, HILIC HPLC method that demonstrates a high degree of efficiency by exploiting these polar interactions between the stationary phase and the target analytes. The chromatography shown in **Figures 7-10** highlight the excellent speed and the high efficiency of this separation while maintaining good retention, peak resolution and peak symmetry.

Peak quantitation was completed by comparison of the peak areas demonstrated in **Figure 1**, to the peak areas of the identified peaks in each of the matrices (**Figures 7-10**). A calibration curve was then made for each analyte to determine the concentration for creatine, creatinine, caffeine, and dicyandiamide, present in each matrix. The amount of each analyte was initially quantitated for a 1 μ L injection of a 2.5 mg/mL solution (**Table 1**). It was then scaled up to be proportional to the amount of one serving size of each of the Pre-Workout powders (**Table 2**).

By employing this HILIC approach, we were able to generate a method that quickly and efficiently separates the key components of the mixture. Using the Luna 3 μ m HILIC 150 x 4.6 mm column we were able to extend a successful separation of chemical standards to analyzing popular market products.

APPLICATIONS

Luna[®] Ordering Information

3 μ m Minibore Columns (mm)				SecurityGuard [™] Cartridges (mm)	
Phases	50 x 2.0	100 x 2.0	150 x 2.0	4 x 2.0*	
HILIC	00B-4449-B0	00D-4449-B0	00F-4449-B0	/10pk	
				AJ0-8328	
				for ID: 2.0-3.0 mm	



3 μ m MidBore [™] and Analytical Columns (mm)					SecurityGuard Cartridges (mm)	
Phases	50 x 3.0	150 x 3.0	100 x 4.6	150 x 4.6	4 x 2.0*	4 x 3.0*
HILIC	00B-4449-Y0	00F-4449-Y0	00D-4449-E0	00F-4449-E0	/10pk	/10pk
					AJ0-8328	AJ0-8329
					for ID: 2.0-3.0 mm	3.2-8.0 mm

5 μ m MidBore Columns (mm)		SecurityGuard Cartridges (mm)	
Phases	150 x 3.0	4 x 2.0*	
HILIC	00F-4450-Y0	/10pk	
		AJ0-8328	
		for ID: 2.0-3.0 mm	

5 μ m Analytical Columns (mm)				SecurityGuard Cartridges (mm)	
Phases	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*	
HILIC	00D-4450-E0	00F-4450-E0	00G-4450-E0	/10pk	
				AJ0-8329	
				for ID: 3.2-8.0 mm	

* SecurityGuard Analytical Cartridges require holder, Part No.: KJ0-4282

guarantee

If Luna analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.



APPLICATIONS

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Luna is a registered trademark, MidBore and SecurityGuard are trademarks of Phenomenex. C4 Cellucor is a registered trademark of Woodbolt Distribution, LLC. N.O. Xplode is a trademark of Gianbia Nutritionals LLC. Bang is a registered trademark of Vital Pharmaceuticals, Inc. Outlift is a trademark of Nutrex Research, Inc. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Dimensions and chromatographic conditions are the same for all columns unless otherwise noted. Comparative separations may not be representative of all applications.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6, 162, 362
CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

© 2016 Phenomenex, Inc. All rights reserved.