


APPLICATION

Lovastatin USP

Overview

The related substances test of the USP monograph outlines the separation of all relevant impurities from Lovastatin. This method was studied and improvements were made to provide higher resolution (Rs) and a faster separation time within allowable adjustments.



Lovastatin

USP Monograph: Lovastatin Details

System Suitability Solution	Dissolve USP Lovastatin RS and USP Lovastatin Related Compound A RS in acetonitrile to obtain a concentration of 2.0 µg/mL of each
Standard Solution	Dissolve USP Lovastatin RS in acetonitrile to obtain a concentration of about 2.0 µg/mL
Test Solution	Dissolve 25 mg of Lovastatin in a 25 mL volumetric flask and dilute to volume with acetonitrile, mix

Column

Size	250 x 4.6 mm
Stationary Phase	5 µm, L7: Octyl silane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod
Temperature	40 °C
Mobile Phase	Acetonitrile and 0.01 M Phosphoric acid (13:7)
Flow Rate	1.5 mL/min
Detection	Spectrophotometer @ 200 nm
Injection	10 µL

Relative Retention with Reference to Lovastatin*

Related Compound A	about 1.3
---------------------------	-----------

System Suitability

Minimum resolution of 6.0 between Lovastatin and Related Compound A

* Retention times, relative retentions, and retardation factors are provided for information only and are not mandatory, no deviation allowance is defined.


Method 1

Original Method within as Described in the USP Monograph

Column: Luna[®] C8(2) 5 µm Fully Porous
Dimensions: 250 x 4.6 mm
Part No.: 00G-4249-E0
Flow Rate: 1.5 mL/min
Sample: 1. Lovastatin
 2. Related Compound A

Elution Time of Last Peak: 10.9 min

Rs Lovastatin and Related Compound A: 12.33



App ID 24118



Method 2


Faster and Higher Resolution Within Allowable Adjustments

Column: Kinetex® Core-Shell C8 5 µm
Dimensions: 250 x 4.6 mm
Part No.: 00G-4608-E0
Flow Rate: 1.5 mL/min
Sample: 1. Lovastatin
 2. Related Compound A

Elution Time of Last Peak: 5.6 min

**Rs Lovastatin
and Related Compound A:** 9.92

Run time reduced by ~5 minutes
with Kinetex Core-Shell Columns!



Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified in Monograph Details Table	As specified
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute change; cannot be reduced to zero	As specified in Monograph Details Table	As specified
Wavelength of Detector	No deviations permitted	200 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	10 µL (as specified)	As specified
Column Temperature	± 10 °C	40 °C (as specified)	As specified
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C8 by C18)	L7 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	250 mm (as specified)	As specified
Column Internal Diameter	Can be adjusted so long as linear velocity if maintained	4.6 mm (as specified)	As specified
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	5 µm (as specified)	As specified
Flow Rate	± 50 % (at given ID)	1.5 mL/min (as specified)	As specified

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25 % to +50 %

Kinetex® Ordering Information

5 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges [‡]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
C8	—	00B-4608-AN	00D-4608-AN	—	AJO-8784 for 2.1 mm ID

5 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
C8	00B-4608-YO	00D-4608-YO	—	AJO-8777 for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
C8	00B-4608-E0	00D-4608-E0	00F-4608-E0	00G-4608-E0	AJO-8770 for 4.6 mm ID

[‡]SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

Luna® Ordering Information

5 µm Microbore and Minibore Columns (mm)								SecurityGuard™ Cartridges (mm)	
Phases	50 x 1.0	150 x 1.0	250 x 1.0	30 x 2.0	50 x 2.0	150 x 2.0	250 x 2.0	4 x 2.0*	
C8(2)	—	00F-4249-A0	—	00A-4249-B0	00B-4249-B0	00F-4249-B0	00G-4249-B0	/10pk	
								AJO-4289	for ID: 2.0-3.0 mm

5 µm MidBore and Analytical Columns (mm)								SecurityGuard™ Cartridges (mm)	
Phases	30 x 3.0	50 x 3.0	150 x 3.0	250 x 3.0	30 x 4.6	50 x 4.6	75 x 4.6	4 x 2.0*	4 x 3.0*
C8(2)	00A-4249-Y0	00B-4249-Y0	00F-4249-Y0	00G-4249-Y0	00A-4249-E0	00B-4249-E0	00C-4249-E0	/10pk	/10pk
								AJO-4289	AJO-4290
								for ID: 2.0-3.0 mm	3.2-8.0 mm

5 µm Analytical and Semi-Prep Columns (mm)					SecurityGuard™ Cartridges (mm)	
Phases	100 x 4.6	150 x 4.6	250 x 4.6	250 x 10	4 x 3.0*	10 x 10 [‡]
C8(2)	00D-4249-E0	00F-4249-E0	00G-4249-E0	00G-4249-N0	/10pk	/3pk
					AJO-4290	AJO-7222
					for ID: 3.2-8.0 mm	9-16 mm

*SecurityGuard™ Analytical Cartridges require holder, Part No.: KJO-4282

[‡]SemiPrep SecurityGuard Cartridges require holder, Part No.: AJO-9281



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to a competing product of the same particle size, similar phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.



APPLICATION

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex and Luna are registered trademarks and MidBore and SecurityGuard are trademarks of Phenomenex.

© 2017 Phenomenex, Inc. All rights reserved.