

# A Comprehensive Screening of Illicit and Pain Management Drugs from Whole Blood Using SPE and LC/MS/MS

## Introduction

Drug analysis from whole blood is gaining popularity due to a more complete measurement of analytes in the biological system. Historically, most drug screens have been performed using immunoassay and other similar technologies. Positive samples are then processed accordingly and analyzed using LC/MS/MS or GC/MS/MS for confirmation and further quantification. Immunoassay methods often lack the specificity and sensitivity necessary for today's complex legal requirements. As a result, LC/MS/MS is a more proficient screening method when an effective pretreatment and sample cleanup method is utilized.

The pretreatment of whole blood samples is an important step in efficiently extracting analytes of interest and maintaining a simplified matrix to analyze. This involves hemolysis to release drugs of interest taken up by the cells as

well as precipitation of proteins. However, determining a single pretreatment option for all classes of drugs can be difficult due to varying chemical properties of the target analytes.

Here we demonstrate a fast, effective technique for the targeted analysis of a broad range of drugs from whole blood utilizing a core-shell HPLC/UHPLC column, SPE, and tandem mass spectrometry. We attempted analysis of thirty nine compounds, comprising opiates (both natural and synthetic), amphetamines, benzodiazepines, phenylpiperidines, a muscle relaxant, and some illicit drugs (PCP and cocaine metabolite). Utilizing advanced sample preparation techniques we simplify a complicated matrix to allow for a fast and successful multi-component analysis by LC/MS/MS.

**Table 1. List of Pain Panel Drugs**

Class	Analyte	Class	Analyte
Benzodiazepines	Alprazolam	Synthetic Opioids	Methadone
	Clonazepam		EDDP
	Diazepam		Fentanyl
	Flunitrazepam		Norfentanyl
	Lorazepam		Meperidine
	Midazolam		Normeperidine
	Nordiazepam		Naloxone
	Oxazepam		Norpropoxyphene
	Temazepam		Propoxyphene
	$\alpha$ -Hydroxyalprazolam		Sufentanil
Alprazolam	Naltrexone		
Opiates	Codeine	Amphetamines	Amphetamine
	Hydrocodone		Methamphetamine
	Hydromorphone		MDMA
	Morphine		MDA
	6-Acetylmorphine (6-MAM)		MDA
	Oxymorphone		
Illicit Drugs	Phencyclidine	Analgesics	Tramadol
	Benzoylcegonine		Carisoprodol
			Buprenorphine
			Norbuprenorphine

## Experimental Conditions

**Table 2. List of Pretreatment Conditions Tested (Not Ranked)**

Option #	Pretreatment Conditions
1	10% Trichloroacetic Acid
2	6% Perchloric Acid
3	Zinc Sulfate + Acetonitrile
4	90:10 Acetonitrile:Methanol
5	50:50 Acetonitrile:Methanol
6	10:90 Acetonitrile:Methanol
7	Zinc sulfate + 90:10 Acetonitrile:Methanol (Most successful. Used for method.)

## SPE Method

<b>Cartridge:</b>	Strata™-X-C 30 mg/3 mL
<b>Part No.</b>	8B-S029-TBJ
<b>Pretreatment:</b>	Add 100 µL 5 % ZnSO <sub>4</sub> ·7H <sub>2</sub> O and lightly vortex. Add 1.5 mL of 90:10 Acetonitrile:Methanol while vortexing. Centrifuge samples at 6000 rpm for 10 min. Add 4 mL of water to the samples to dilute the amount of organic prior to SPE
<b>Condition:</b>	1 mL Methanol
<b>Equilibrate:</b>	1 mL Water
<b>Wash 1:</b>	1 mL 0.1 % Formic Acid
<b>Wash 2:</b>	1 mL 30 % Methanol
<b>Dry:</b>	3 to 4 min at high vacuum (~10" of Hg)
<b>Elute:</b>	2 x 500 µL (2 aliquots of 500 µL) Ethyl acetate:Isopropanol:Ammonium Hydroxide (70:20:10)
<b>Blow Down:</b>	To dryness under Nitrogen at 40-45 °C
<b>Reconstitute:</b>	With 500 µL of 85:15 (A:B) of LC mobile phase

## LC Conditions

<b>Column:</b>	Kinetex® 2.6 µm Biphenyl		
<b>Dimensions:</b>	50 x 3 mm		
<b>Part No.:</b>	00B-4622-Y0		
<b>Mobile Phase:</b>	0.1% Formic Acid in Water / 0.1% Formic Acid in Methanol		
<b>Gradient:</b>	<b>Time</b>	<b>A</b>	<b>B</b>
	0.0	90	10
	2.5	0	100
	3.5	0	100
	3.51	90	10
	5	90	10
<b>Flow Rate:</b>	0.7 mL/min		
<b>Injection Volume</b>	10 µL		
<b>Temperature:</b>	Ambient		

### MS/MS Parameters


Detection: AB Sciex 4000 QTRAP®, positive polarity

ESI CAD Gas = High; GS1=65, GS2 = 50; curtain gas = 10



ESI spray voltage = 4000 V; Temp = 550°C

## Results

**Figure 1. Representative Chromatogram of Drug Panel**




**Figure 2. Selectivity of Isobaric Pair Morphine and Hydromorphone**



**Hydromorphone**

$C_{17}H_{19}NO_3$

MW: 285.34




**Morphine**

$C_{17}H_{19}NO_3$

MW: 285.34

### Figure 3. Comparison of Select Pretreatment Extracts Prior to SPE




A higher percent of acetonitrile results in a clearer extract

A= 10:90 Acetonitrile:Methanol  
B= 50:50 Acetonitrile:Methanol  
C= 90:10 Acetonitrile:Methanol




## Figure 4. Comparison of the effects of various pretreatment options on amphetamine

Chromatograms are overlaid with time shift to provide clarity.



## Figure 5. Comparison of the effects of various pretreatment options on benzoylecgonine

Chromatograms are overlaid with time shift to provide clarity.



## Summary of Sample Pretreatment Procedure

Different classes of compounds responded better with different pretreatment options


- Zinc Sulfate + Acetonitrile combination was most effective at extracting the benzodiazepines, synthetic opioids, amphetamines, analgesics, and other illicit drugs
- 90:10 Acetonitrile:Methanol produced the best recoveries for the natural opiates

Based on these findings, a pretreatment combination of Zinc Sulfate and 90:10 Acetonitrile:Methanol was determined to be the best option to extract the majority of analytes of interest

### Table 3. Calibration Curve Parameters

Class	Analytes	Corr. Coeff.	Class	Analytes	Corr. Coeff.	
Benzodiazepines	Alprazolam	0.9968	Opiates	Codeine	0.9919	
	Clonazepam	0.9967		Oxycodone	0.9932	
	Diazepam	0.9960		Hydromorphone	0.9942	
	Flunitrazepam	0.9959		Morphine	0.9949	
	Lorazepam	0.9948		6-Acetylmorphine	0.9948	
	Midazolam	0.9974		Hydrocodone	0.9979	
	Amphetamines	Nordiazepam	0.9968	Synthetic Opioids	Methadone	0.9911
		Oxazepam	0.9966		EDDP	0.9960
		Temazepam	0.9957		Fentanyl	0.9934
$\alpha$ -Hydroxylalprazolam		0.9953	Norfentanyl		0.9935	
Amphetamine		0.9981	Meperidine		0.9906	
Methamphetamine	0.9912	Normeperidine	0.9954			
MDMA	0.9912	Naloxone	0.9923			
MDA	0.9963	Norpropoxyphene	0.9944			
MDEA	0.9923	Propoxyphene	0.9933			
Analgesics	Tramadol	0.9930	Sufentanil		0.9972	
	Carisoprodol	0.9931	Naltrexone		0.9987	
	Buprenorphine	0.9910	Others		Phencyclidine	0.9906
Norbuprenorphine	0.9943	Benzoylcegonine			0.9902	

**Figure 6. Calibration Curve for Various Analytes Over 200 Fold Concentration (5-1000 ng/mL)**



**Table 4. Method Precision and Accuracy Data Based on Replicate Quality Control Samples**

Analyte	Class	Expected Conc. ng/mL (Low)	%RSD (Low)	% Accuracy (Low)	Expected Conc. ng/mL (High)	%RSD (High)	% Accuracy (High)
Alprazolam	Benzodiazepines	20	10	108	200	12	104
Clonazepam		20	9	114	200	11	107
Diazepam		20	10	97	200	12	103
Flunitrazepam		20	7	112	200	7	105
Lorazepam		20	15	108	200	10	111
Midazolam		20	7	115	200	4	88
Nordiazepam		20	11	101	200	13	103
Oxazepam		20	6	108	200	12	105
Temazepam		20	7	105	200	9	99
$\alpha$ -Hydroxyalprazolam		20	6	88	200	11	91
Codeine	Opiates	20	10	92	200	9	87
Oxycodone		20	4	95	200	2	93
Hydromorphone		20	6	85	200	14	97
Hydrocodone		20	7	105	200	9	99
Morphine		20	8	91	200	10	86
Methadone	Synthetic Opioids	20	10	110	200	5	105
EDDP		20	10	98	200	2	94
6-MAM		20	7	100	200	7	100
Fentanyl		20	9	115	200	5	90
Norfentanyl		20	12	95	200	4	100
Meperidine		20	7	105	200	7	103
Normeperidine		20	9	103	200	10	102
Naloxone		20	7	118	200	3	111
Norpropoxyphene		20	9	100	200	14	90
Propoxyphene		20	12	111	200	5	101
Sufentanil		20	8	98	200	7	89
Naltrexone		20	4	113	200	11	108
Amphetamine	Amphetamines	20	9	107	200	11	107
Methamphetamine		20	10	115	200	3	96
MDMA		20	13	111	200	8	92
MDA		20	8	102	200	7	101
MDEA		20	16	107	200	3	105
Tramadol	Analgesics	20	4	105	200	3	96
Carisoprodol		20	8	106	200	9	100
Buprenorphine		20	12	104	200	11	101
Norbuprenorphine		20	6	105	200	13	106
Phencyclidine	Others	20	7	110	200	4	92
Benzoylcegonine	Others	20	10	104	200	5	101

## Discussion

A combination of zinc sulfate with acetonitrile and methanol worked for the majority of analytes including opiates and synthetic opioids, benzodiazepines, and amphetamines.

The Kinetex 2.6  $\mu\text{m}$  Biphenyl 50 x 3 mm column provided quick separation of all drugs of abuse analytes in a 5 minute cycle time with a total elution window of ~2 minutes (**Figure 1**).

- Excellent resolution and good selectivity of the opiate isobaric species was obtained, such as morphine and hydromorphone (**Figure 2**).
- Fast elution time as a result of a ballistic gradient led to the potential of even greater throughput with the use of multiplex or dual-stream systems.

Strata-X-C is effective for all the classes of pain panel drugs that are typically monitored in whole blood testing (amphetamines, benzodiazepines, opiates, and other drugs of abuse).

- Successfully removes interfering endogenous compounds such as proteins, which improves sensitivity and also extends column lifetime.
- Allows for sample concentration which greatly increases the sensitivity of assays and allows for older, less sensitive mass spectrometry systems to be effectively used for drugs of abuse testing.

Good linearity ( $R^2 \geq 0.990$ ) was demonstrated for all drugs of abuse compounds over a wide concentration range (**Table 3** and **Figure 6**).

Precision and accuracy were within suitable limits ( $\pm 15\%$ ) (**Table 4**).

## Conclusion

We have developed an effective pretreatment and SPE cleanup method for whole blood followed by targeted LC/MS/MS analysis. Zinc sulfate with an acetonitrile and methanol combination provided the best response for the majority of analytes tested, and the Strata-X-C SPE sorbent chemistry worked well at removing endogenous matrix interferences while retaining a wide range of compounds. The pretreatment/SPE method showed robustness over a dynamic range of

concentrations, demonstrating good precision and accuracy. The sample cleanup procedure described here can greatly extend column lifetime and reduce maintenance of the MS system. The use of the Kinetex 2.6  $\mu\text{m}$  Biphenyl provides the necessary selectivity and resolving power to analyze all target drug molecules in a total cycle time of about 5 minutes, including baseline separation of critical isobaric ions.



## References

1. S Huq, S Sadjadi, and S Countryman, "Quantitative Bio Analysis of the Most Commonly Used Pain Medications in Urine Using a Reliable Sample Preparation Technique in Combination With an API 5000 LC-MS-MS", Mass Spec Application for Clinical Laboratory Conference, 2013
2. Dalsgaard et al, Quantitative Analysis of 30 Drugs in Whole Blood by SPE and UHPLC-TOF-MS. (2013) *J of Forensic Science and Criminology* 1(1):101

### Trademarks

Strata-X is a trademark and Kinetex is a registered trademark of Phenomenex. API 4000 is a trademark and QTRAP is a registered trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

© 2014 Phenomenex, Inc. All rights reserved.