

TN-1120

APPLICATIONS

LC/MS/MS Analysis of Chloramphenicol from Various Matrices using Strata[™]-X SPE and a Kinetex[®] 2.6 μm C18 HPLC/UHPLC Column


Several food matrices, including shrimp, milk, and honey were analyzed by LC/MS/MS for chloramphenicol (CAP).

This method has been optimized for the following AB SCIEX mass spectrometers: API 3200[™], 3200 QTRAP[®], API 4000[™], 4000 QTRAP, Triple Quad[™] 4500, QTRAP 4500, Triple Quad 5500, and QTRAP 5500.

Analytes

Chloramphenicol: C₁₁H₁₂Cl₂N₂O₅

CAS: 56-75-7 MW = 311.015 amu



For AB SCIEX mass spectrometer users, this method can be instantly implemented by installing the iMethod[™] Test. iMethod Tests are verified across several laboratories and contain everything you will need to start running samples including sample preparation recommendations, consumables, LC running conditions, optimized MRM parameters, reporting templates, and complete method documentation.

Visit www.phenomenex.com/iMethod for more information about available iMethod Tests.

Materials and Methods

Extracting Chloramphenicol from Shrimp

Sample Preparation

1. Homogenize ~100 g of thawed shrimp using a blender or tissue homogenizer.
Note: A small amount of water may need to be added to facilitate homogenization of the tissue.
2. Weigh out 5 g of homogenized shrimp and transfer to a 15 mL polypropylene tube.
3. Add 50 μL of d₅-Chloramphenicol internal standard (I.S.) solution.
4. Mix thoroughly using a vortex mixer to ensure adequate distribution of the I.S. throughout the homogenate.
5. Add 2 mL of water.
6. Mix well using a vortex mixer.
7. Add 5 mL of ethyl acetate.
8. Transfer the tubes to a mechanical shaker and shake rigorously for 30 minutes.
9. Centrifuge at 7000 rpm for 15 minutes.
10. Transfer the supernatant to a new 15 mL polypropylene test tube, reserving the pelleted shrimp homogenate tissue.
11. Add 5 mL of ethyl acetate to the pelleted shrimp homogenate tissue from step 10 and repeat the extraction process.
12. Combine the resulting ethyl acetate supernatant with the extraction previously obtained in step 10.
13. Dry the combined supernatant extracts under nitrogen gas at 55 °C.
14. Reconstitute with 300 μL of methanol, and then dilute to 10 mL with water.

Note: The sample is ready for solid phase extraction.

TN-1120

APPLICATIONS

Solid Phase Extraction (SPE)

SPE can be done using a Phenomenex Strata® SPE Vacuum Manifold*

- 12-Position Manifold (Part No. AH0-6023)
- 24-Position Manifold (Part No. AH0-6024)

* Manifold includes: vacuum-tight glass chamber, polypropylene lid with gasket, bleed valve and gauge, stopcock valves, collection racks, polypropylene needles.

Sorbent:	Strata™-X, 60 mg/3 mL
Part No.:	8B-S100-UBJ
Condition:	3 mL Methanol
Equilibrate:	3 mL Water
Load:	10 mL pre-treated sample
Wash:	3 mL Water
Dry:	5 min at 10 in. of Hg
Elute:	3 mL Ethyl acetate
Dry Down:	Dry down @ 55 °C under a stream of nitrogen
Reconstitute:	Reconstitute in 0.5 mL of water/acetonitrile (80:20)

The sample can be filtered through a Phenex™ 0.45 µm filter and then transferred to a Verex™ autosampler vial for LC/MS/MS analysis.

Ethyl acetate should be used as a rinse solution to prevent chloramphenicol contamination (for example, used to rinse needle from blow-down station or interior of the stopcock that is positioned into the manifold).

Comments on the Extraction and SPE of Chloramphenicol

Chloramphenicol has a tendency to exhibit carry-over effects. Extreme caution and good laboratory techniques should be used in order to minimize these effects. In particular, make sure that the SPE manifold and blow-down apparatus are thoroughly cleaned before and after every use (methanol is an appropriate solvent for cleaning these devices). In addition, it is always a good practice to run an extracted solvent blank (processed exactly as the normal samples, but lacking the shrimp tissue) if contamination of any of the solvents is suspected.

TN-1120

APPLICATIONS

Extraction Techniques for Chloramphenicol in other Foodstuffs

Sample Preparation (Honey)

Based on a method by Nieland et al.¹

1. Dissolve 5 g of honey in 20 mL of water.
2. Add 100 µL of internal standard (D₅-CAP 1 µg/mL).
3. Mix using a vortex mixer until uniform suspension.
4. Extract with 20 mL ethyl acetate in an ultrasonic bath for 5 minutes.
5. Centrifuge for 10 minutes at 4000 rpm.
6. Transfer the upper phase (ethyl acetate) into a glass test tube.
7. Evaporate until dry under a stream of nitrogen and gentle heating (55 °C).
8. Dissolve the residue in 1 mL of water.
9. Filter using a Phenex™ 0.45µm filter and then transfer to a Verex™ autosampler vial for LC/MS/MS analysis.

Sample Preparation (Royal Jelly)

1. Measure 2 g of royal jelly.
2. Add 100 µL of internal standard (D₅-CAP 1 µg/mL).
3. Add 20 g of Septra™ Bulk Florisil® (Part No. 04G-4411).
4. Homogenize with a glass rod.
5. Load homogenized mix onto a Strata® FL-PR (Florisil) 5 g/20 mL Tube (Part No. 8B-S013-LEG)
6. Flush vial, glass rod and column with 250 mL of ethyl acetate/cyclohexane (50:50).
7. Evaporate to dryness under a stream of nitrogen and gentle heating (55 °C).
8. Dissolve the residue in 1 mL of water.
9. Filter using a Phenex 0.45µm filter and then transfer to a Verex autosampler vial for LC/MS/MS analysis.

Sample Preparation (Fish and Seafood)

Based on a method by Effkemann et al.²

1. Measure and homogenize 5 g of fish or seafood tissue.
2. Add 100 µL of internal standard (D₅-CAP 1 µg/mL).
3. Add 10 mL of ethyl acetate.
4. Homogenize for 1 minute using an ULTRA-TURRAX®.
5. Centrifuge for 10 minutes at 3000 rpm.
6. Transfer the upper phase (ethyl acetate) into a glass test tube.
7. Evaporate until dry under a stream of nitrogen and gentle heating (55 °C).
8. Dissolve the residue in 1 mL of water.
9. Add 1 mL of hexane and mix.
10. Centrifuge for 10 minutes at 3000 rpm.
11. Remove the lower phase (water) and then filter it through a Phenex 0.45µm filter and then transfer to a Verex autosampler vial for LC/MS/MS analysis.

TN-1120

APPLICATIONS

Sample Preparation (Milk)

Based on a method by Hormazabal et al.³

1. Measure 2 mL of milk sample.
2. Add 100 μ L of internal standard (D₅-CAP 1 μ g/mL).
3. Extract with 5 mL of acetonitrile.
4. Centrifuge for 10 minutes at 4000 rpm.
5. Wash supernatant with 5 mL of chloroform.
6. Centrifuge for 10 minutes at 4000 rpm.
7. Discharge upper phase (water).
8. Dry lower phase under a stream of nitrogen.
9. Dissolve residue in 100 μ L of methanol.
10. Add 3 mL of water and 0.5 mL disodium hydrogen phosphate (0.5 M, pH 6.0).
11. Condition a Strata[®] SDB-L 200 mg/6 mL SPE cartridge (Part No. 8B-S014-FCH) with 2 mL methanol followed by 2 mL water.
12. Load sample onto SPE cartridge.
13. Wash twice with 2 mL of water and 0.3 mL of water/methanol (50:50).
14. Elute twice with 0.3 mL of methanol.
15. Evaporate until dry under a stream of nitrogen and gentle heating (55 °C).
16. Reconstitute in 1 mL of water.
17. Filter using a Phenex[™] 0.45 μ m filter and then transfer to a Verex[™] autosampler vial for LC/MS/MS analysis.

Sample Preparation (Milk Powder)

Based on a method by Hormazabal et al.³


1. Measure 5 g of milk powder.
2. Add 100 μ L of internal standard (D₅-CAP 1 μ g/mL).
3. Dilute with 25 mL of water.
4. Homogenize for one minute using an ULTRA-TURRAX[®].
5. Continue with step 3 of Sample Preparation (Milk).

TN-1120

APPLICATIONS

HPLC Methodology

Figure 1. Extracted ion chromatograms for the lowest level of chloramphenicol extracted from spiked shrimp tissue (0.05 µg CAP/Kg of shrimp tissue)



The signal-to-noise ratio for the quantifier ion (321.2/152) was 531.8:1. The signal-to-noise ratio for the qualifier ion (321.2/257) was 59.9:1.

LC/MS/MS Conditions

Column: Kinetex® 2.6 µm C18
Dimensions: 50 x 2.1 mm
Part No.: 00B-4462-AN
Mobile Phase: A: Water
B: Acetonitrile
Gradient:

Time (min)	B (%)
0.00	5
2.01	95
4.00	95

Flow Rate: 0.4 mL/min
Injection Volume: 25 µL
Temperature: 25 °C
Detection: MS/MS, ESI negative (ESI-)
Sample: 1. Chloramphenicol
2. Chloramphenicol-d₅

TN-1120

APPLICATIONS

MS/MS Detection	
Optimized for API 4000™ and 4000 QTRAP® LC/MS/MS systems	
TurboV™ source	
Negative polarity	

Source/Gas Parameters	
CUR:	10 psi
IS:	-4500 V
TEM:	450 °C
GS1:	45 psi
GS2:	50 psi
Ihe:	On
CAD:	5

Compound MRM Parameters and Retention Times	
Dwell Time:	150 ms for analytes and internal standards
MRM Pause Time:	5 ms
Q1 Resolution:	UNIT
Q3 Resolution:	UNIT

Retention times were collected using the Agilent® 1200 LC system.

API 3200™ and 3200 QTRAP® system parameters

Analyte	Q1	Q3	DP	EP	CE	CXP	Agilent RT (min)	Shimadzu® RT (min)
Chloramphenicol 1	321.0	151.9	-30	-4	-24	-4	2.09	2.08
Chloramphenicol 2	321.2	257.2	-30	-4	-14	-6	2.09	2.08
D ₅ -Chloramphenicol	326.1	157.0	-30	-4	-26	-5	2.08	2.07

API 4000™, 4000 QTRAP, Triple Quad™ 4500, and QTRAP 4500 system parameters

Analyte	Q1	Q3	DP	EP	CE	CXP	Agilent RT (min)	Shimadzu® RT (min)
Chloramphenicol 1	321.0	152.0	-50	-10	-24	-4	2.09	2.08
Chloramphenicol 2	321.2	257.2	-50	-10	-16	-6	2.09	2.08
D ₅ -Chloramphenicol	326.1	157.0	-50	-10	-26	-5	2.08	2.07

Triple Quad 5500 and QTRAP 5500 system parameters

Analyte	Q1	Q3	DP	EP	CE	CXP	Agilent RT (min)	Shimadzu® RT (min)
Chloramphenicol 1	321.0	151.9	-50	-4	-24	-4	2.09	2.08
Chloramphenicol 2	321.0	257.0	-50	-4	-14	-6	2.09	2.08
D ₅ -Chloramphenicol	326.1	157.0	-50	-4	-26	-5	2.08	2.07

Results Table for the Standards

Analyte	Expected Concentration (ng/mL)	Number of Values	Mean	Standard Deviation	% CV	% Accuracy
CAP	0.500	2	0.485	0.027	5.6	97.0
CAP	5.000	4	5.878	0.113	1.9	117.6
CAP	25.000	4	23.075	0.814	3.5	92.3
CAP	50.000	4	43.096	0.944	2.2	86.2
CAP	200.000	4	215.431	5.793	2.7	107.7

References

1. B. Nieland et al.: 'Detection of a range of antibiotics in honey by Liquid Liquid Extraction followed by LC/MS/MS' poster presented at ASMS Conference on Mass Spectrometry (2004) Nashville, Tennessee, USA
2. S. Effkemann: 'Bestimmung von Chloramphenicol in Fischen mittels RP-HPLC mit massenspektrometrischer Detektion' Prüfvorschrift LC-MS002 LAVES (2005) Cuxhaven, Germany
3. V. Hormazabal, M. Yndestad: J. Liq. Chrom. & Rel. Technol. 24 (2001) 2477

TN-1120

APPLICATIONS

Ordering Information



For AB SCIEX mass spectrometer users, this method can be instantly implemented by installing the iMethod™ Test. iMethod Tests are verified across several laboratories and contain everything you will need to start running samples including sample preparation recommendations, consumables, LC running conditions, optimized MRM parameters, reporting templates, and complete method documentation.

Visit www.phenomenex.com/iMethod for more information about available iMethod Tests.

Complete iMethod™ Kit*

Description	Part No.
iMethod Test for Chloramphenicol	KHO-8985

* Kit includes Kinetex® 2.6 µm C18 HPLC/UHPLC column, SecurityGuard™ ULTRA cartridges and holder, in-line filter, Phenex™ Syringe Filters, Verex™ Vial kit, Sure-Lok™ UHPLC Nut and tightening tool, Sure-Lok Fingertight Nut and wrenches, and Strata™-X SPE cartridges.

Strata-X SPE

Sorbent Mass	Part No.	Unit
Tube		
30 mg	8B-S100-TAK	1 mL (100/box)
30 mg	8B-S100-TBJ	3 mL (50/box)
60 mg	8B-S100-UBJ	3 mL (50/box)
100 mg	8B-S100-EBJ	3 mL (50/box)
100 mg	8B-S100-ECH	6 mL (30/box)
200 mg	8B-S100-FBJ	3 mL (50/box)
200 mg	8B-S100-FCH	6 mL (30/box)
500 mg	8B-S100-HBJ	3 mL (50/box)
500 mg	8B-S100-HCH	6 mL (30/box)
Giga™ Tube		
500 mg	8B-S100-HDG	12 mL (20/box)
1 g	8B-S100-JDG	12 mL (20/box)
1g	8B-S100-JEG	20 mL (20/box)
2 g	8B-S100-KEG	20 mL (20/box)
5 g	8B-S100-LFF	60 mL (16/box)

Additional sizes available. Contact your Phenomenex Sample Preparation Specialist for additional information.

Strata® SDB-L SPE

Sorbent Mass	Part No.	Unit/Box
Tube		
100 mg	8B-S014-EAK	1 mL (100/Box)
200 mg	8B-S014-FBJ	3 mL (50/Box)
200 mg	8B-S014-FCH	6 mL (30/Box)
500 mg	8B-S014-HBJ	3 mL (50/Box)
500 mg	8B-S014-HBJ	6 mL (30/Box)
1g	8B-S014-JCH	6 mL (30/Box)
Giga Tube		
10 g	8B-S014-MFF	60 mL (16/Box)

Strata Florisil® SPE

(pesticide residue grade)

Sorbent Mass	Part No.	Unit/Box
Tube		
500 mg	8B-S013-HBJ	3 mL (50/Box)
500 mg	8B-S013-HCH	6 mL (30/Box)
1g	8B-S013-JCH	6 mL (30/Box)
Giga Tube		
1 g	8B-S013-JEG	20 mL (20/Box)
2 g	8B-S013-KDG	12 mL (20/Box)
5 g	8B-S013-LEG	20 mL (20/Box)
10 g	8B-S013-MFF	60 mL (16/Box)

Kinetex Core-Shell HPLC/UHPLC Columns

1.7 µm Minibore Columns (mm)

Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
C18	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJ0-8782 for 2.1 mm ID

SecurityGuard™
ULTRA Cartridges†

2.6 µm Minibore Columns (mm)

Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
C18	00A-4462-AN	00B-4462-AN	00D-4462-AN	00F-4462-AN	AJ0-8782 for 2.1 mm ID

SecurityGuard™
ULTRA Cartridges†

2.6 µm MidBore™ Columns (mm)

Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
C18	00A-4462-Y0	00B-4462-Y0	00C-4462-Y0	00D-4462-Y0	00F-4462-Y0	AJ0-8775 for 3.0 mm ID

SecurityGuard™
ULTRA Cartridges†

2.6 µm Analytical Columns (mm)

Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
C18	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	AJ0-8768 for 4.6 mm ID

SecurityGuard™
ULTRA Cartridges†

†SecurityGuard ULTRA cartridges require holder, Part No. AJ0-9000



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and comparable dimensions, return the product with comparative data within 45 days for a FULL REFUND.

TN-1120 APPLICATIONS

Australia

t: 02-9428-6444
f: 02-9428-6445
auinfo@phenomenex.com

Austria

t: 01-319-1301
f: 01-319-1300
anfrage@phenomenex.com

Belgium

t: 02 503 4015 (French)
t: 02 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: (800) 543-3681
f: (310) 328-7768
info@phenomenex.com

Denmark

t: 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: 09 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: 01 30 09 21 10
f: 01 30 09 21 11
franceinfo@phenomenex.com

Germany

t: 06021-58830-0
f: 06021-58830-11
anfrage@phenomenex.com

India

t: 040-3012 2400
f: 040-3012 2411
indiainfo@phenomenex.com

Ireland

t: 01 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: 051 6327511
f: 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: 030-2418700
f: 030-2383749
nlinfo@phenomenex.com

New Zealand

t: 09-4780951
f: 09-4780952
nzinfo@phenomenex.com

Norway

t: 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: (800) 541-HPLC
f: (310) 328-7768
info@phenomenex.com

Sweden

t: 08 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: 01625-501367
f: 01625-501796
ukinfo@phenomenex.com

United States

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com

All other countries: Corporate Office USA

t: (310) 212-0555
f: (310) 328-7768
info@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex and Strata are registered trademarks of Phenomenex. Strata-X, Phenex, Giga, MidBore, SecurityGuard, Verex, Sepra and Sure-Lok are trademarks of Phenomenex. iMethod, Cliquid, and QTRAP are registered trademarks of AB SCIEX, Ltd. Triple Quad API 3200, API 4000 and TurboV are trademarks of AB SCIEX, Ltd. Shimadzu is a registered trademark of Shimadzu Scientific Instruments. Florisil is a registered trademark of U.S. Silica Company. ULTRA-TURRAX is a registered trademark of IKA.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

© 2012 Phenomenex, Inc. All rights reserved.