BEFORE USING COLUMN

PLEASE READ IMPORTANT TIPS

Starting with Aeris WIDEPORE

Because of the reduced hydrophobicity compared to fully porous 300Å columns, one should start gradients with reduced organic concentrations to increase retention and improve peak shape of polar proteins and peptides. Shallower gradients compared to fully porous columns may also be appropriate.

Starting with Aeris PEPTIDE

For increased resolving power, use a longer column (preferably a 250 mm) to maximize separation. For specific peptide map applications, one can also use different gradient slope segments to “stretch” or “compress” regions of a peptide map.

Flow Rate and Column Length

Start method development using column dimensions and flow rates similar to existing protein or peptide separation methods. For higher throughput, use shorter column lengths and higher flow rates for minimal impact on resolution. For maximum resolution, use longer columns and shallower gradients.

Column ID	Flow Rate
4.6 mm ID Columns	0.8-2.0 mL/min
2.1 mm ID Columns	0.2-0.5 mL/min

HPLC System Optimization

Aeris 2.6 µm and 3.6 µm core-shell columns operate comfortably within the pressure limits of conventional HPLC systems and meet or exceed the performance of sub-2 µm fully porous particle columns on UHPLC systems. To maximize the benefits of your Aeris Core-Shell columns, one should investigate:

- Minimizing peak dispersion by reducing the system dwell volume between the injector and detector.
- Optimizing detector settings by adjusting the scan rate and/or time constant to the fastest practical setting such that signal-to-noise ratio (S/N) is not adversely affected.

The Core-Shell Performance Enhancement Kit is a full hardware set designed to minimize system dwell volume.

Questions? Ask AerisBallGirl@phenomenex.com

EASILY OPTIMIZE YOUR HPLC SYSTEM TO MAXIMIZE THE PERFORMANCE OF AERIS CORE-SHELL COLUMNS

- Increase method efficiency, resolution, and detection
- Make connections quickly and reliably – every time, on any HPLC system
- Improve results for sensitive and demanding applications

The fittings and tubing used in this kit are carefully chosen to minimize dwell volume and reduce band broadening. Combined with your Aeris Core-Shell column, this kit will provide reliable connections and quality performance every time.

View the Installation Video at

www.phenomenex.com/OptimizeAeris


Ordering Information

Core-Shell Performance Enhancement Kit

Part No.	Description	Unit
AQ0-8892	Core-Shell Performance Enhancement Kit, Includes: PEEKsil™ Tubing, Fittings and Tool	ea

INSTALLATION DIAGRAM


Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Aeris and SecurityGuard are trademarks of Phenomenex. PEEKsil is a trademark of SGE.

© 2014 Phenomenex, Inc. All rights reserved.

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© 2014 Phenomenex, Inc. Tutti i diritti riservati.