

APPLICATIONS

A Comparison of Various Kinetex[®] C18 Phases

Zeshan Aqeel and Jeff Layne
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Introduction

When performing reversed phase LC separations, there is a natural tendency of method developers to first test a C18 phase, the most commonly used phase in LC applications. Today, there are many different C18 phases available to method developers so choosing a simple C18 column may prove to be more challenging than one had originally thought.


This technical note explores the different C18 phases within the Kinetex core-shell family of HPLC/UHPLC columns to determine the differences as well as to provide a general guideline of how to select the most appropriate C18 phase for a variety of applications.

Materials and Methods

To provide a streamlined method development approach, the following cases were studied to determine the most appropriate C18 phase for each application:


- Acid, Base, Neutral Panel
- Highly Polar Analytes
- Highly Polar Bases (Catecholamines)
- Highly Polar Bases (Nucleotides)
- Hydrophobic Neutrals and Acids (Cannabinoids)

Three different Kinetex core-shell C18 phases were screened under the same conditions as well as optimized conditions in some cases:




Kinetex Polar C18

Combined C18 and polar modified surface that provide polar and non-polar retention alongside 100 % aqueous stability



Kinetex C18

Balanced C18 phase that provides the highest degree of hydrophobic selectivity relative to the other Kinetex phases



Kinetex XB-C18

Unique C18 phase that yields increased hydrogen bonding with hydrophobic selectivity, resulting in improved peak shape for basic compounds and increased retention of acidic compounds



Results and Discussion


I. Acid, Base, Neutral Panel using a Formic Acid Gradient

- In general, gradients will reduce selectivity differences between phases as the gradient rate itself will play a large role in analyte elution
- But, particularly when using formic acid and other weak buffers, differences in ion-exchange and polar interactions will still be apparent, giving subtle (or not so subtle) differences in selectivity
- These differences will be less apparent when using more strong buffer salts like phosphate, or when using ion-pairing agents (e.g. hexane sulfonic acid) or ion-masking agents (TEA; TFA)


LC Column: As specified
Dimensions: 50 x 4.6 mm
Mobile Phase: A: 0.1 % Formic acid in Water
 B: 0.1 % Formic acid in Acetonitrile
Gradient: Time (min) B (%)
 0 5
 0.5 5
 10 95
 11 95
Flow Rate: 1.85 mL/min
Injection: 1 µL
Temperature: 30 °C

Sample: 1. Pyridine 8. Phenol
 2. Acetaminophen 9. Triprolidine
 3. Sulfathiazole 10. Prednisolone
 4. Quinidine 11. Nortriptyline
 5. Quinidine Impurity 12. 5-Methylsalicylaldehyde
 6. Acetubotolol 13. Hexanophenone
 7. Chlorpheniramine


Kinetex® 2.6 µm Polar C18



Kinetex 2.6 µm C18



Kinetex 2.6 µm XB-C18




II. Highly Polar Analytes – Water-Soluble Vitamins


- Take advantage of the aqueous stability of Kinetex® Polar C18 to perform your analyses in 100 % aqueous conditions to maximize retention of polar analytes
- Traditional C18 phases like Kinetex C18 or XB-C18 should not be used below 2-3 % organic to maintain phase stability

LC Column: As specified
Dimensions: 50 x 4.6 mm
Mobile Phase: 20 mM Potassium phosphate, pH 3.5 / Methanol (97:3) (except where noted)
Flow Rate: 1 mL/min
Injection: 1 µL
Temperature: 35 °C
Detection: UV @ 210 nm
Sample: 1. Thiamine
 2. Nicotinamide
 3. Pyridoxal
 4. Pyridoxine


Kinetex 2.6 µm Polar C18 (mobile phase = 100% 20 mM Potassium phosphate, pH 3.5)




Kinetex 2.6 µm Polar C18



Kinetex 2.6 µm C18



Kinetex 2.6 µm XB-C18




III. Highly Polar Bases – Catecholamines

- Take advantage of the aqueous stability on Kinetex[®] Polar C18 to perform your analyses in 100 % aqueous conditions to maximize retention of polar analytes
- Traditional C18 phases like Kinetex C18 or XB-C18 should not be used below 2-3 % organic to maintain phase stability


LC Column: As specified
Dimensions: 50 x 4.6 mm
Mobile Phase: 20 mM Ammonium acetate, pH 6.9 / Methanol (98:2) (except where noted)
Flow Rate: 1 mL/min
Injection: 1 µL
Temperature: Ambient
Detection: UV @ 280 nm

Sample: 1. Norepinephrine
 2. Epinephrine
 3. Dopamine
 4. Serotonin
 5. 5-Hydroxy-3-Indoleacetic Acid
 6. 4-Hydroxy-3-Methoxyphenylacetic Acid


Kinetex 2.6µm Polar C18 (mobile phase = 100 % 20 mM Ammonium acetate, pH 6.9)




Kinetex 2.6µm Polar C18



Kinetex 2.6µm C18



Kinetex 2.6µm XB-C18




IV. Highly Polar Bases – Nucleotides

- Take advantage of the aqueous stability on Kinetex[®] Polar C18 to perform your analyses in 100 % aqueous conditions to maximize retention of polar analytes. The trade-off can be an increase in tailing, as is apparent below, but we include the data as a logical extreme; one would expect to optimize their method to find the right balance of retention and peak shape
- The key is that operating in highly aqueous conditions is only possible with the Kinetex Polar C18 phase


LC Column: As specified
Dimensions: 50 x 4.6 mm
Mobile Phase: 20 mM Ammonium acetate, pH 6.9 / Methanol (97:3) (except where noted)
Flow Rate: 1 mL/min
Injection: 1 µL
Temperature: Ambient
Detection: UV @ 280 nm

Sample: 1. Thiourea
 2. 5-Fluorocytosine
 3. Adenine
 4. Thymine
 5. Guanosine-5-Monophosphate


Kinetex 2.6µm Polar C18 (mobile phase = 100 % 20 mM Ammonium acetate, pH 6.9)




Kinetex 2.6µm Polar C18



Kinetex 2.6µm C18



Kinetex 2.6µm XB-C18




V. Hydrophobic Neutrals and Acids – Cannabinoids

- Even very hydrophobic molecules will display selectivity differences on different C18 phases due to structural (e.g. positional isomers) or polar functional groups
- In this case, although selectivity was similar overall, the traditional Kinetex[®] C18 provided better selectivity for one critical pair (CBN and CBGA), which is interesting because they are totally different structurally one would think could easily separate


LC Column: As specified
Dimensions: 50 x 4.6 mm
Mobile Phase: A: 0.1 % Formic acid in Water
 B: 0.1 % Formic acid in Acetonitrile
Gradient: Time (min) B (%)
 0 60
 16 100
Flow Rate: 0.5 mL/min
Injection: 1 µL
Temperature: 50 °C
Detection: UV @ 230 nm

Sample: 1. CBDV 7. CBN
 2. CBDVA 8. CBGA
 3. THC-V 9. D9-THC
 4. CBD 10. D8-THC
 5. CBG 11. CBC
 6. CBDA 12. THCA-A


Kinetex 2.6 µm C18



Kinetex 2.6 µm Polar C18



Kinetex 2.6 µm XB-C18



Conclusions

- You will find differences in selectivity between the Kinetex[®] C18, Kinetex XB-C18, and Kinetex Polar C18 stationary phases based upon the differences in surface chemistry
- Under identical running conditions, the differences may be subtle – these are after all still all C18 phases – but subtle differences can make the difference between resolution and co-elution, or lead to reversals in elution order that can help to resolve a target analyte from matrix interferences
- With the newest Kinetex Polar C18 phase, however, you can operate in 100 % aqueous conditions, whereas the traditional C18 phases should not go below 2-3 % organic. This opens the

potential to generate significant changes in retention for polar molecules.

- A thorough method development project should involve screening each of these phases, but if only one column could be screened, it is recommended that it be the Kinetex Polar C18, with its ability to operate in 100 % aqueous conditions, followed by either the Kinetex C18 or Kinetex XB-C18 to complement it
- The Kinetex Polar C18, with its aqueous stability and special surface modification, may provide a benefit for analyzing polar analytes, whereas the standard Kinetex C18 or XB-C18 phases, may be advantageous for separations based upon primarily hydrophobic differences
- Although it was not included in this study, our Kinetex EVO C18 is the obvious choice for operating in alkaline mobile phases, or when you need improved peak shape for tailing basic analytes, especially when using weakly-buffered mobile phases like 0.1 % formic acid

Kinetex Ordering Information

1.3 µm Minibore Columns (mm)		
Phases	30 x 2.1	50 x 2.1
C18	00A-4515-AN	00B-4515-AN

1.7 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges [†]
Phases	30 x 3.0	50 x 3.0	100 x 3.0	3/pk
XB-C18	00A-4498-YO	00B-4498-YO	00D-4498-YO	AJO-8775
C18	—	00B-4475-YO	00D-4475-YO	AJO-8775

for 3.0 mm ID

1.7 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
EVO C18	00B-4726-AO	00D-4726-YO	00F-4726-AO

1.7 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	—	00B-4726-AN	00D-4726-AN	00F-4726-AN	AJO-9298
XB-C18	00A-4498-AN	00B-4498-AN	00D-4498-AN	00F-4498-AN	AJO-8782
C18	00A-4475-AN	00B-4475-AN	00D-4475-AN	00F-4475-AN	AJO-8782

for 2.1 mm ID

2.6 µm Analytical Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	3/pk
EVO C18	—	00B-4725-EO	—	00D-4725-EO	00F-4725-EO	AJO-9296
Polar C18	—	00B-4759-EO	—	00D-4759-EO	00F-4759-EO	AJO-9530
XB-C18	—	00B-4496-EO	00C-4496-EO	00D-4496-EO	00F-4496-EO	AJO-8768
C18	00A-4462-EO	00B-4462-EO	00C-4462-EO	00D-4462-EO	00F-4462-EO	AJO-8768

for 4.6 mm ID

2.6 µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges [†]
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	—	00B-4725-YO	—	00D-4725-YO	00F-4725-YO	AJO-9297
Polar C18	—	00B-4759-YO	—	00D-4759-YO	00F-4759-YO	AJO-9531
XB-C18	00A-4496-YO	00B-4496-YO	00C-4496-YO	00D-4496-YO	00F-4496-YO	AJO-8775
C18	00A-4462-YO	00B-4462-YO	00C-4462-YO	00D-4462-YO	00F-4462-YO	AJO-8775

for 3.0 mm ID

2.6 µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges [†]
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4725-AN	00B-4725-AN	—	00D-4725-AN	00F-4725-AN	AJO-9298
Polar C18	00A-4759-AN	00B-4759-AN	—	00D-4759-AN	00F-4759-AN	AJO-9532
XB-C18	00A-4496-AN	00B-4496-AN	00C-4496-AN	00D-4496-AN	00F-4496-AN	AJO-8782
C18	00A-4462-AN	00B-4462-AN	00C-4462-AN	00D-4462-AN	00F-4462-AN	AJO-8782

for 2.1 mm ID

More Ordering Information →

[†] SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

APPLICATIONS

Kinetex[®] Ordering Information

2.6 Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
XB-C18	00B-4496-A0	00F-4496-E0	00F-4496-A0

3.5 µm Analytical Columns (mm)			SecurityGuard ULTRA Cartridges [‡]
Phases	100 x 4.6	150 x 4.6	3/pk
XB-C18	00B-4496-A0	00D-4796-E0	AJO-8768

for 4.6 mm ID

5 µm Minibore Columns (mm)					SecurityGuard [™] ULTRA Cartridges [‡]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJO-9298
XB-C18	00A-4605-AN	00B-4605-AN	00D-4605-AN	—	AJO-8782
C18	00A-4601-AN	00B-4601-AN	00D-4601-AN	00F-4601-AN	AJO-8782

for 2.1 mm ID

5 µm MidBore [™] Columns (mm)				SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	00B-4633-Y0	00D-4633-Y0	00F-4633-Y0	AJO-9297
XB-C18	00B-4605-Y0	00D-4605-Y0	00F-4605-Y0	AJO-8775
C18	00B-4601-Y0	00D-4601-Y0	00F-4601-Y0	AJO-8775

for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJO-9296
XB-C18	00B-4605-E0	00D-4605-E0	00F-4605-E0	00G-4605-E0	AJO-8768
C18	00B-4601-E0	00D-4601-E0	00F-4601-E0	00G-4601-E0	AJO-8768

for 4.6 mm ID

5 µm Semi-Preparative Columns (mm)			SecurityGuard SemiPrep Cartridges ^{***}
Phases	150 x 10	250 x 10	10 x 10
EVO C18	00F-4633-N0	00G-4633-N0	AJO-9306
C18	00F-4601-N0	00G-4601-N0	AJO-9278

for 10 mm ID

5 µm Axia [™] Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges [‡]
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	15 x 21.2
EVO C18	00B-4633-P0-AX	00D-4633-P0-AX	00F-4633-P0-AX	00G-4633-P0-AX	AJO-9304
XB-C18	00B-4605-P0-AX	00D-4605-P0-AX	00F-4605-P0-AX	00G-4605-P0-AX	AJO-9145
C18	00B-4601-P0-AX	00D-4601-P0-AX	00F-4601-P0-AX	00G-4601-P0-AX	AJO-9145

for 21.2 mm ID

5 µm Axia Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges [‡]
Phases	50 x 30	100 x 30	150 x 30	250 x 30	15 x 30
EVO C18	00B-4633-U0-AX	00D-4633-U0-AX	00F-4633-U0-AX	00G-4633-U0-AX	AJO-9305
XB-C18	00B-4605-U0-AX	00D-4605-U0-AX	00F-4605-U0-AX	00G-4605-U0-AX	AJO-9204
C18	00B-4601-U0-AX	00D-4601-U0-AX	00F-4601-U0-AX	00G-4601-U0-AX	AJO-9204

for 30 mm ID

‡ SecurityGuard ULTRA Cartridges require holder, Part No.: AJO-9000

* PREP SecurityGuard Cartridges require holder, Part No.: AJO-8223

** PREP SecurityGuard Cartridges require holder, Part No.: AJO-8277

*** SemiPrep SecurityGuard Cartridges require holder, Part No.: AJO-9281

Australia
t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China
t: +86 400-606-8099
f: +86 (0)22 2532-1033
phen@agela.com

Denmark
t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India
t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico
t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico
t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain
t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

**All other countries
Corporate Office USA** 
t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



If Kinetex core-shell analytical columns do not provide at least an equivalent separation as compared to a competing column of the same particle size, similar phase and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

Trademarks

Kinetex is a registered trademark and SecurityGuard, Axia, and MidBore are trademarks of Phenomenex.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7,674,383.

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2017 Phenomenex, Inc. All rights reserved.

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com