


APPLICATIONS

Robust Separation of Hydroxychloroquine and Chloroquine in Hydroxychloroquine Sulfate Tablets Using the Kinetex® 5 µm C18

Zeshan Aqeel, J Preston, Ryan Splitstone, and Phil Koerner
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Monograph: Hydroxychloroquine Sulfate
pK_a: 9.67
LogP: 3.87

Overview

In this application is the comparison of two particle morphologies of similar L1 (C18) phase in the separation of Hydroxychloroquine and related impurity Chloroquine under the USP conditions for assay of hydroxychloroquine sulfate tablets. The application demonstrates the potential method improvements that can be achieved per the allowable adjustments outlined in USP General Chapter <621> relative to the original column and conditions referenced in the monograph.

USP Monograph: Hydroxychloroquine Sulfate Tablet Assay

Standard Stock Solution	Dissolve 1.0 mg/mL of USP Hydroxychloroquine Sulfate RS in <i>Diluent</i>
Standard Solution	0.05 mg/mL of USP Hydroxychloroquine Sulfate RS from <i>Standard Stock solution</i> in <i>Mobile Phase</i>
System Suitability Stock Solution	Dissolve 1.0 mg/mL of Chloroquine Phosphate in Methanol
System Suitability Solution	Transfer 5.0 mL of the <i>System Suitability Stock solution</i> into 100-mL volumetric flask, add 5.0 mL of <i>Standard Stock solution</i> , and dilute to volume with <i>Mobile Phase</i>
Diluent	Mixture of Methanol:Water (1:1)

Column


Size	Method 1: 250 x 4.6 mm, Method 2: 250 x 4.6 mm
Stationary Phase	Method 1: Fully Porous 5 µm C18, Method 2: Kinetex 5 µm C18
Temperature	26 °C
Mobile Phase (premixed)	A = Methanol, Acetonitrile, water, and phosphoric acid (100:100:800:2). Mix and add 96 mg of sodium 1-pentanesulfonate in the resulting solution, and filter Mobile Phase
Isocratic	Isocratic: (100:100:800:2, Methanol: Acetonitrile: Water: Phosphoric Acid) Total Run Time: 20 min
Flow Rate	1.0 mL/min
Detector	UV @ 254 nm
Injection Volume	10 µL of System Suitability solution and Standard solution

System Suitability – System Suitability solution and Standard solution

Sample: Standard solution and System Suitability solution:

- Resolution (Rs): NLT 1.8 between Chloroquine and Hydroxychloroquine for System Suitability solution
- Relative Standard Deviation: NMT 1.5 % for Standard solution (5 replicate injections)

USP Tailing Factor, <621> Chromatography**



Symmetry factor (A_S): Also known as the “tailing factor”, of a peak is calculated by:

$$A_S = W_{0.05}/2f$$


where W_{0.05} is the width of the peak at 5% height and f is the distance from the peak maximum to the leading edge of the peak, the distance being measured at a point 5% of the peak height from the baseline.

** General Chapter <621> “Chromatography” in United States Pharmacopeia 40 National Formulary 35 (USP 40-NF 35, United States Pharmacopeial Convention, Rockville, Maryland, 2017), p. 6.

APPLICATIONS

Method 1


Standard solution: 0.05 mg/mL Hydroxychloroquine Sulfate



App ID: 25859

Method 2


Standard solution: 0.05 mg/mL Hydroxychloroquine Sulfate



App ID: 25857

Method 1

System Suitability solution: 0.05 mg/mL Hydroxychloroquine Sulfate & 0.05 mg/mL Chloroquine




App ID: 258560

APPLICATIONS

Method 2

mAU System Suitability solution: 0.05 mL/mg Hydroxychloroquine Sulfate & 0.05 mL/mg Chloroquine



Column: Kinetex[®] 5 µm C18 Core-Shell
Dimensions: 250 x 4.6 mm
Part No.: 00G-4601-E0
Flow Rate: 1.0 mL/min
Sample: 1. Chloroquine
 2. Hydroxychloroquine

Average Res: 6.16
Average Res RSD: 0.28 %

Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified	As specified
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute adjustment; cannot be reduced to zero	As specified	As specified
Wavelength of Detector	No deviations permitted	254 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	20 µL (as specified)	10 µL (Allowed)
Column Temperature	± 10 °C	Ambient (as specified)	26 °C (Allowed)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	250 mm (as specified)	250 mm (Allowed)
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.6 mm (as specified)	4.6 mm (Allowed)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	5 µm (as specified)	5 µm (Allowed)
Flow Rate	± 50 % (at given ID)	1.0 mL/min (as specified)	1.0 mL/min**

*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25 % to +50 %.

**Maintained volumetric flow rate at the given column ID.

Allowable Column Adjustments: L/dp Ratio -25 % to 50 %

Column	Length (mm)	ID (mm)	dp (µm)	L/dp	Allowable Range
Fully Porous	250	4.6	5	50,000	37,500 - 75,000
Core-Shell	250	4.6	5	50,000	ALLOWED

Method Summary and Comparison

	Method 1	Method 2
Column	Fully Porous 5 µm C18	Kinetex 5 µm C18
Hydroxychloroquine Average Rt	9.8 min	6.2 min
Hydroxychloroquine Average Rs	3.42	6.16
Standard Solution Hydroxychloroquine Peak Area	391.6	412.8
Standard Solution Hydroxychloroquine Peak Area RSD (n=5)	2.26 %	0.28 %
Backpressure (Bar)	176	190

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/ChatNow.

Australia

t: +61 (0)2-9428-6444
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 (12) 881 0121
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☎ **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark and BE-HAPPY is a trademark of Phenomenex.

Disclaimer

Chromatograms may not be representative of all separations.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.