



presston™ 100

YOUR NEWEST SOLUTION TO INCREASED PRODUCTIVITY



POSITIVE PRESSURE MANIFOLD
for Streamlined Sample Preparation

 **phenomenex**[®]
...breaking with traditionSM



www.phenomenex.com/Presston

Do More with Presston™ 100

Process 96-well plates and 1 mL, 3 mL, or 6 mL tubes using a single manifold*

Sleek, Low Profile Design

Width: 13", Depth: 13.25", Height: 17" (open), 13.5" (closed)

Safely Process Samples

Protect operators from solvents

Never Lose Pressure

Always maintain a tight seal between the manifold and tubes/96-well plate

Determine Your Operating Pressure

Easily read your operating pressure with the easy to view pressure gauge

Easily Load Samples

Moveable locator plate makes sample loading and cleaning easy



Phenomenex warrants that for a period of 12 months following delivery, the Presston 100 Positive Pressure Manifold you have purchased will perform in accordance with the published specifications and will be free from defects in materials or workmanship. In the event that the Presston 100 Positive Pressure Manifold does not meet this warranty, Phenomenex will repair or replace defective parts.

Please visit www.phenomenex.com/Presston for complete warranty information.

*Tube adapter kit required

Ordering Information




Part No.	Description	Unit Price
AH0-9334	Presston 100 Positive Pressure Manifold, 96-Well Plate	
AH0-9342	Presston 100 Positive Pressure Manifold, 1 mL Tube Complete Assembly	
AH0-9347	Presston 100 Positive Pressure Manifold, 3 mL Tube Complete Assembly	
AH0-9343	Presston 100 Positive Pressure Manifold, 6 mL Tube Complete Assembly	

The Presston 100 96-Well Positive Pressure Manifold can also process 1, 3, and 6 mL tubes using the following adapter kits

Presston 100 Tube Adapter Kits (for AH0-9334)

Part No.	Description	Unit Price
AH0-9344	1 mL Tube Adapter Kit	
AH0-9345	3 mL Tube Adapter Kit	
AH0-9346	6 mL Tube Adapter Kit	

Presston 100 Pairs Well with the Following Premium Sample Preparation Solutions



Visit www.phenomenex.com/Presston to learn more



Your newest solution to increased productivity

Positive Pressure Manifold for Streamlined Sample Preparation

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auiinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Strata is a registered trademark and Novum, Phree, Impact, β -Gone, and Presston are trademarks of Phenomenex.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

© 2016 Phenomenex, Inc. All rights reserved.